



Iridex® LIO Plus

An Easy to Use Single-Mirror Laser Indirect
Ophthalmoscope That Is Ergonomic and
Provides **Excellent Peripheral Visualization**,
Treatment Flexibility, Consistency, and Reliability



iridex.com/lio

The gold standard of LIOs

For many years, the Iridex® LIO Plus, a single-mirror laser indirect ophthalmoscope, has been considered the **gold standard of LIOs**. Retina specialists and teaching institutions around the world use it for its excellent visualization, ease of use, and reliability.



.....

- The Iridex LIO Plus is indicated for transpupillary retinal photocoagulation. It facilitates observation and photocoagulation of peripheral retinal pathologies, including diabetic retinopathy and retinopathy of prematurity.
- Provides excellent peripheral visualization and consistent laser uptake as far as the ora serrata
- It is ergonomic, intuitive, and ideal to use in patients who are best examined and treated in a supine position
- Has a lightweight headset with size adjustments to optimize fit and comfort
- Permits independent positioning of the laser spot within the illuminated field, or simultaneous adjustment of both the laser spot and illuminated field within the user's visual field
- The long laser depth of focus tolerates a wide range of working distances
- The 810 nm Large Spot model can be used with the OcuLight SLx Laser to deliver transpupillary thermotherapy for the treatment of intraocular tumors.^{1,2}
- The single-wavelength models are compatible with 532 nm, 577 nm, and 810 nm Iridex lasers
- A dual-wavelength model can be used with both 532 nm and 810 nm Iridex lasers for added utility
- The halogen illumination source offers consistent illumination and great color rendering
- Its integrated laser filters provide protection for the physician during laser use
- Includes a travel case for safe storage and ease of transportation to satellite offices



“Treating comfortably, with minimal ergonomic burden, is an important way to help consistency of results across treatments. The LIO Plus, with its intuitive design and practical ensure consistency of results across treatments. The LIO Plus, with its intuitive design and practical features, plays a large role in my laser treatment protocol – It is proving to be effective for my most challenging peripheral cases, and my colleagues and I have been pleased with the results.”

Shawn Kavoussi, MD • Texas Retina Center, Houston, United States



Compatibility, Specifications, and Ordering Information

Product Number	Wavelength	Spot Size*	Laser Compatibility	Connector
30903 - H500	532 nm / 810 nm	360 µm	OcuLight GL, GLx, TX, SL, SLx, IQ 532, DioVet	Resistive
13152 - H500	810 nm	360 µm	OcuLight SL, SLx, DioVet	Resistive
13153 - H500 (Large Spot)	810 nm	1400 µm	OcuLight SL, SLx	Resistive
65515 - H500	532 nm	360 µm	IQ 532	RFID
65900 - H500	577 nm	360 µm	IQ 577	RFID

Dimensions

Width	7.0 cm / 2.75 in
Length	13.5 cm / 5.31 in
Height	12.3 cm / 4.84 in
Weight	239 g / 0.52 lbs

Ophthalmoscope

Headband circumference adjustment range	520 mm to 640 mm
Illumination field sizes**	50 mm, 40 mm, 20 mm
Interpupillary adjustment	48 mm - 74 mm
Ophthalmoscope filters	Clear, red-free, amber

To request a quote or a demo, speak with your local sales representative or email us at retina@iridex.com.

* May vary with refractive power. Retinal spot when treating an emmetropic patient eye with 20 D lens.

** At aerial image plane.

¹Hasanreisoglu M, Saktanasate J, Schwendeman R, Shields JA, Shields CL. Indocyanine green-enhanced transpupillary thermotherapy for retinoblastoma: Analysis of 42 tumors J Pediatr Ophthalmol Strabismus, 2015;52(6):348-54.

²Shields CL, Santos MC, Diniz W, Gunduz K, Mercado G, Cater JR, Shields JA. Thermotherapy for retinoblastoma Arch Ophthalmol, 1999;117(7):885-93.



EC REP Emergo Europe
Prinsessegracht 20,
2514 AP The Hague, The Netherlands



Contact Iridex today to learn more



Iridex • 1212 Terra Bella Avenue • Mountain View, CA 94043 • U.S.A.
800.388.4747 (U.S.A. only) • 1.650.940.4700 • retina@iridex.com • iridex.com/lto.de



© 2022 Iridex. All rights reserved. Iridex, the Iridex logo, OcuLight, IQ 577, IQ 532 and trademarks or registered trademarks of Iridex.
88290-EN Rev. A 05.2022